

Upgrading to Phys-X Advanced (ver. 7.0)

If you are upgrading from an older version of Phys-X to version 7.0 please follow these instructions to import your existing Prescriptions, Custom Exercises & Protocols.

1. Install Phys-X 7.0 to its own directory. (If you are unsure, simply use the default when installing.)
2. Run the new version of the program. From the main screen select the "**Setup...**" menu and choose "**Import Data from Previous Version...**"
3. In the dialog box that appears click the **Browse** button.
4. Navigate to the older Phys-X application directory & select the ***.cfg** file. (example: Physx6.cfg) Click the **Open** button.
5. If checked, the **Favor Photos** option will set your imported Protocols to using photographs, if available for the exercise. (We recommend you leave it checked.)
6. Once the previous *.cfg file has been selected, click the **Import Data** button.
7. The progress of importing will be displayed. If you have a large number of prescriptions and custom exercises this process may take several minutes.
8. After you verify that all your files have been successfully converted, you can un-install the older version.